

SUMMIT MAX



MAXimum Strength

MAXimum Comfort

MAXimum Durability

MAXimum Safety

500 lb. Capacity

Extra Wide Seat

Reinforced Footrest

 **SUMMIT LIFTS, INC.**
Taking you to new heights

866-378-6648

SUMMIT **WAY**

23" wide, high-backed contoured seat helps ensure a comfortable, secure ride.

Features:

Industry leading weight capacity
Spacious comfortable contoured seat
Smooth starts and stops
Operates on household current

Specifications:

Capacity-500 lb.
Top seat height- 20"
Bottom seat height- 23"
Travel-Up to 20'
Track-Extruded Aluminum
Motor-1/2 hp 90 VDC
Electric Supply-110 VAC grounded outlet
Cable-Aircraft cable rated at 4,200 lb. breaking strength
Finish-Durable off white powder coat finish,
choice of 3 upholsteries
Controls-Constant pressure pendant control, and
choice of 2 hand or foot activated call controls
(24 volt) or wall-mounted call-send controls
Warranty-2 year warranty on all component parts
and drive train

Safety:

- Constant-pressure controls stop the lift the instant you release the control
- Limit switches automatically shut unit off at the top and bottom of the stairs
- Seat swivels 180 degrees toward top landing for safe and easy entry and exit
- Safety switch prevents lift from running when the seat is not in proper position

The seat swivels at the top landing making it easier to get on and off the lift. Safety sensors prevent the lift from running when the seat is not in the proper position.

Heavy-duty armrests help provide more secure transfers.

Reinforced footrest for added stability.



Due to overall size of lift, a minimum 36" wide staricase is recommended.



Choose from three colors of upholstery:



SUMMIT LIFTS, INC.
Taking you to new heights

18505 E. 163rd St., Lake Winnebago, MO 64034
Toll Free: 866-378-6648
E-mail: info@summitlifts.com
Visit us at www.summitlifts.com