

# Stairlifts



STAIRLIFTS / WHEELCHAIR LIFTS / ELEVATORS / RAMPS / PATIENT LIFTS / BATHROOMS



# EHL S



*Accessibility Solutions Since 1991*

# If you have difficulty getting up and down your stairs, you're not alone.

Every year thousands of people look for a solution to the problems associated with difficulty navigating the stairs in their homes. The decision to install a stairlift is sometimes a difficult one, but when compared with the alternatives many people find a stairlift is the easiest and best option.

## Stairlift Basics



Stairlifts easily fold when not in use creating space on stairs. They are battery-operated and charged on a standard outlet.



Rails are mounted to the treads of the stairs and don't require a wall for support.



The offset swivel seat mechanism provides a safe, comfortable exit/entry at the top of the stairs. As a safety feature the lift will not operate in this position.



The two wireless remotes allow caregiver operation, two user operation, and parking the chair up or down so chair can be stored out of view.



Safety features are built into the stairlift. The footrest sensors safely stop the stairlift when the slightest obstruction is detected.



**Whether your staircase is Straight or Curved — EHL has a wide selection of stairlifts that fit the most elaborate staircases.**

**We deliver solutions that fit your specific needs.**

**Our stairlifts are battery-operated, easily installed, reliable, and safe.**


- **Straight and Curved Staircases**
- **Manual or Automatic Folding Rail**
- **Power Swivel Seat**
- **Upholstery Choices**
- **Painted Track**
- **Larger Chair and Footrest**



Straight



Curved



Curved Stair Configuration

## Straight Stairlifts



### ELAN—Bruno

- 300 lbs. capacity
- Continuous charging rail
- Installs within 5" of wall
- Power/manual folding rail
- Power options for swivel and footrest available



### ELITE—Bruno

- 400 lbs. capacity
- Continuous charging rail
- Power/manual folding rail
- Options/Upgrades available
- Outdoor model available
- Power options for swivel and footrest available



### PINNACLE—Harmar

- 350 lbs. capacity
- Continuous charging rail
- 600 lbs. heavy duty model option



### SIMPLICITY (+) Handicare

- 310 lbs. capacity
- Continuous charging rail
- Power/manual folding rail



### 1000—Handicare

- 350 lbs. capacity
- Power/manual folding rail
- Options/Upgrades available

## Curved Stairlifts



### ELITE CRE—Bruno

- 400 lbs. capacity
- Factory custom built for your staircase
- Options/Upgrades available



### HW 10—Hawle

- 350 lbs. capacity
- Dual rail tubular system
- Outdoor model available

## Visit the EHLS Showroom

Test-Ride a Stairlift in our operational showroom located at:

**210 W Campus Drive  
Arlington Heights, IL**

*"I love the stair lifts and am glad I can get around my home again. The installer and sales people were especially nice and informative."*

*— Regina W, Chicago*

## Key Stairlift Options:



Up/Down Rocker Switch Controls Swivel Seat

(Not available with the Larger Seat Option)



Power Footrest



Larger Footrest



Larger Chair

Park Positions on Curved Rails  
90/180 Degrees  
Provides Ease & Safety



Manual or Power Folding Rails



# Why EHLS

## Local and Experienced:

- Family-owned and locally operated since 1991
- We serve the entire Chicago Metro Area, including SE Wisconsin & NW Indiana
- We install the equipment we sell and our technicians and installers are factory trained

## Accessible Remodeling Solutions:

- Your one stop shop for all accessibility solutions including: stairlifts, elevators, wheelchair lifts, and bathroom/home remodeling.

## Service:

- We stand behind our work and will be there if your equipment needs repair
- We provide same or next day service in most cases

## Accessibility Showroom:

- We offer a completely operational showroom, open to the public to help customers test equipment and ensure that it is right for them
- The showroom is equipped with Straight and Curved Stair Lifts, a Wheelchair lift, a Home Elevator, a Ramp, Transfer/Ceiling lifts, a Roll-In Shower and an Accessible Bathroom

## We Make House Calls:

- We will come to your home and assess the options that best meet YOUR needs
- Our Home Assessments are free

## We Stock Lifts/Equipment:

- Our large warehouse allows us to keep a selection of lifts and equipment in stock for quick turn-around
- Our leading equipment suppliers are U.S. and North American based
- Our primary equipment manufacturers of stairlifts, wheelchair lifts and elevators are based in the Midwest

### Stairlift Snapshot:

STRAIGHT STAIRLIFTS	Folding Rail		Power Swivel Seat	Larger Footrest	Weight Capacity	Upholstery Choices	Painted Track
	Manual	Auto					
BRUNO ELAN	✓	✓	✓	✓	300 lbs.		✓
BRUNO ELITE	✓	✓	✓	✓	400 lbs.	✓	✓
HARMAR PINNACLE		✓		✓	350 lbs.	✓	✓
HANDICARE SIMPLICITY (+)	✓		✓	✓	310 lbs.	✓	
HANDICARE 1000	✓		✓	✓	350 lbs.	✓	
CURVED STAIRLIFTS							
BRUNO ELITE CRE			✓	✓	400 lbs.	✓	✓
HAWLE HW10			✓	✓	350 lbs.	✓	✓



Serving the greater Chicagoland area, Northwest Indiana, and Southeast Wisconsin.

**Call for a FREE on-site assessment or planning/design consultation.**

**Extended Home Living Services, Inc.**

Visit our Accessibility Showroom

210 W Campus Drive • Suite B • Arlington Heights IL • 60004

Direct (847) 215-9490 | (630) 717-4445 | (773) 775-6122

Fax (847) 590-1728

www.ehls.com



Division of EHLS