

Transfer Lift Systems



STAIRLIFTS / WHEELCHAIR LIFTS / ELEVATORS / RAMPS / PATIENT LIFTS / BATHROOMS



EHLS[®]



Accessibility Solutions Since 1991

Transfer Lift Systems

Ceiling Mounted Lifts (CMLs) are ideal solutions to assist in lifting and transferring those with mobility challenges.



CMLs can:

- Provide a safe and secure way to transfer – without risk
- Protect caregivers from strain when lifting and transferring
- Allow a single caregiver to manage the transfer safely
- Provide more stability than a rolling patient lift
- Accommodate even the most limited living spaces

Permanent Lifts

- Continuous track system that travels within or between rooms
- Battery operated and charged while on the track
- Powered travel allows the caregiver to lift the person up and move along the track

Portable Lifts

- Easily moveable between the tracks and used wherever transfer is needed
- Ideal for when continuous track installation is not an option



Visit the EHLS Showroom

Test-Ride a CML in our operational showroom located at:

**210 W Campus Drive
Arlington Heights, IL**

Rail Systems


With a ceiling mounted lift or overhead lift system, you can lift and transfer with minimal effort and maximum safety.

Our staff will help to design the right solution for your needs.

Contact us for an in-home assessment with one of our accessibility specialists.

"My only regret is not having gotten our portable ceiling lifts sooner. We waited too long and my wrists, shoulders and back could take no more. I highly recommend this product for ALS patients; it is the only safe way to transfer someone who no longer has use of the arms and legs. I beg of you to get the word out so others get the relief they need. Everyone involved were most supportive. Our sincere appreciation for all your help in getting the job done as soon as possible."

— Nancy and Rich W, Hoffman Estates



The lift unit is only part of a well-functioning overhead lift system.

A complete and flexible rail system is just as important. We offer a variety of options for a range of functions and installation types including:

- **Straight Rails**
Ability to move all around the house with added turntables
- **Complex Rails**
Utilizing curves and carpentry to facilitate door access
- **X-Y Gantry Rails**
A three piece system with two parallel fixed rails and perpendicular moving rail that makes for easy adjustment and mobility



Slings

Sling styles vary based on the needs and abilities of the user. Slings are more encompassing when a user needs more support and are less able to assist the caregiver.



With over 20 years of experience, an EHLS Accessibility Specialist can help you make the right decision. Call today for your FREE in-home consultation!

Why EHLS®

Since 1991,

Extended Home Living Services has been a leader in the design and installation of accessibility lifts in homes and public buildings. We have installed hundreds of lifts in houses, schools, churches, restaurants, and offices. Whether it's new construction or a renovation, EHLS has a complete line of equipment that can be integrated with your design to meet your accessibility needs and ADA requirements. EHLS would be happy to review plans or schedule an on-site visit.

Local and Experienced:

- Family-owned and locally operated since 1991
- We serve the entire Chicago Metro Area, including SE Wisconsin & NW Indiana
- We install the equipment we sell and our technicians and installers are factory trained



Accessibility Showroom:

- We offer a completely operational showroom, open to the public to help customers test equipment and ensure that it is right for them
- The showroom is equipped with Straight and Curved Stair Lifts, a Wheelchair lift, a Home Elevator, a Ramp, Transfer/Ceiling lifts, a Roll-In Shower and an Accessible Bathroom

We Make House Calls:

- We will come to your home and assess the options that best meet YOUR needs
- Our Home Assessments are free

We Stock Lifts/Equipment:

- Our large warehouse allows us to keep a selection of lifts and equipment in stock for quick turn-around
- Our leading equipment suppliers are U.S. and North American based
- Our primary equipment manufacturers of stairlifts, wheelchair lifts and elevators are based in the Midwest

Accessible Remodeling Solutions:

- Your one stop shop for all accessibility solutions including stairlifts, elevators, wheelchair lifts, and bathroom/home remodeling.

Service:

- We stand behind our work and will be there if your equipment needs repair
- We provide same or next day service in most cases



Serving the greater Chicagoland area, Northwest Indiana, and Southeast Wisconsin.

Call for a FREE on-site assessment or planning/design consultation.

Extended Home Living Services, Inc.

Visit our Accessibility Showroom

210 W Campus Drive • Suite B • Arlington Heights IL • 60004

Direct (847) 215-9490 | (630) 717-4445 | (773) 775-6122

Fax (847) 590-1728

www.ehls.com



Division of EHLS